



Streamlined HPLC Method Development

A step-by-step protocol for fast method optimization

INTRODUCTION

This guide outlines a structured, step-by-step method development process grounded in selectivity screening, scientific decision-making and proven tools, including ACE® Method Development Kits from Avantor®. Whether you're analyzing small molecules, polar analytes, biomolecules or complex mixtures, this systematic strategy allows you to reduce development time, increase success rates and create high-quality separations that hold up to real-world variability.

STEP 1: Define the analytical objective

Before you start to select columns, solvents or settings, take a moment to clearly define your objective and what success will look like. This information will influence all subsequent decisions.

Ask yourself:

- What analytes must be separated?
- What is the sample matrix? (e.g., biological, environmental, pharmaceutical)
- What level of resolution, sensitivity or quantification is required?
- What is the final application? (e.g., research, regulatory filing, QC testing)
- What detection method will be used? (e.g., UV, MS, ELSD)

When you know the critical pairs, detection limits and method purpose, you will be better able to determine the most appropriate separation mode and screening strategy.

STEP 2: Choose the best column screening kit for your needs

Column chemistry can be the biggest influence on selectivity in HPLC. Rather than sticking with one stationary phase and optimizing around it, successful method development can also start by testing different chemistries side by side.

Avantor ACE Method Development Kits allow scientists to evaluate multiple phase chemistries in parallel, under identical conditions. This allows you to isolate the effects of stationary phase selectivity and make rapid progress toward achieving workable separations.

Here are recommended kits, based on analyte type:

Analyte type	Recommended kit	Includes
Aromatics, drugs, isomers	ACE Advanced Kit	C18, C18-AR, C18-PFP
Polar compounds, pH-sensitive	ACE Extended Kit	SuperC18, C18-Amide, CN-ES
Highly polar, ionic species	ACE HILIC Kit	HILIC-A, HILIC-B, HILIC-N
Basic and neutral analytes, fast methods	ACE UltraCore Kit	UltraCore SuperC18, UltraCore SuperPhenylHexyl
Peptides, proteins	ACE Bioanalytical 300 Å Kit	C18-300, C4-300, Phenyl-300

Each kit contains two or three columns with distinct interaction mechanisms. This will help maximize your chances of separating difficult analytes without needing to alter solvents, gradients or instrumentation.

STEP 3: Establish initial screening conditions

You should start your screening with consistent, universal conditions so that any differences in separations can be attributed solely to different stationary phases.

Here is an example setup for reversed-phase or general screening:

- Mobile phase A: water with 0.1% formic acid
- Mobile phase B: acetonitrile with 0.1% formic acid
- Gradient: 5–95% B over 10 min
- Flow rate: 0.5 ml/min
- Temperature: 30–40 °C
- Detection: 210 or 254 nm, or MS where applicable

By keeping everything except columns consistent, you can more accurately evaluate the impact of column chemistry on analyte retention, selectivity and peak shape without concern that it was influenced by other factors.

STEP 4: Run side-by-side column screens

After you have chosen which columns you want to evaluate, run the same sample on each column in your kit using the same conditions. Pay attention to:

- Retention order
- Baseline separation of critical pairs
- Peak shape (tailing, fronting)
- Co-elutions or unresolved zones
- Total run time and peak distribution

Even subtle changes in bonded phase chemistry can create dramatic improvements in resolution. For example:

- **ACE C18-AR** improves separation of aromatic compounds through π - π interactions.
- **ACE C18-PFP** enhances separation of halogenated or electron-rich molecules through dipole interactions.
- **ACE CN-ES** separates polar compounds with improved retention of dipolar species.
- **ACE HILIC-N** enables separation of very polar, water-soluble analytes under highly organic conditions.

This screening process can help you identify which column helps you achieve your desired selectivity with minimal need for downstream adjustments.

STEP 5: Optimize around the best column

After you have identified the best column for your needs, you can begin to fine-tune your method. Here are some adjustments to consider:

- Gradient profile: Adjust slope and range to improve peak spacing.
- Buffer pH: Alter analyte ionization and retention.
- Organic modifier: Methanol vs. acetonitrile can shift selectivity.
- Temperature: Higher temperatures reduce viscosity and may sharpen peaks.
- Flow rate: Can affect resolution and analysis time.

For the greatest accuracy, make only one change at a time. Be sure to document all changes to provide method traceability and testing.

STEP 6: Confirm robustness and reproducibility

Robust methods hold up under normal variability. Before you lock your method, take these steps and compare the results to make sure they are consistent:

- Repeat with fresh columns from different lots.
- Vary flow rate $\pm 10\%$, temperature ± 5 °C or pH ± 0.2 units.
- Run on different days and analysts.
- Test in relevant matrices (e.g., plasma, tablets, water).

Building robust methods helps ensure your work is suitable for transfer, scale-up or QC deployment.

WHY USE AVANTOR ACE METHOD DEVELOPMENT KITS?

ACE Method Development Kits from Avantor are made to meet three principles that matter to scientists:

- **Selectivity by design.** Each kit contains complementary phases selected to expose differences in retention mechanisms. This helps accelerate your path to determining baseline separation.
- **Ultra-inert silica.** All ACE columns use ultra-inert base silica and proprietary bonding. This helps reduce silanol activity and delivers superior peak shape, even for basic compounds.
- **Expert support.** Avantor provides free method development assistance from experienced chromatographers. Whether you're unsure which kit to choose or need help optimizing pH or gradient conditions, our expert help is available to you.

AN EXAMPLE: SEPARATING STRUCTURALLY RELATED COMPOUNDS

A scientist working on a formulation containing structurally similar analgesics found poor resolution on a traditional C18 column. Screening with an ACE Advanced Kit revealed:

- ACE C18: Strong retention but co-elution of key peaks
- ACE C18-PFP: Slight selectivity improvement
- ACE C18-AR: Full resolution with excellent peak shape and minimal tailing

The selected column was then optimized for gradient and pH to reduce development time by more than 50% compared to previous efforts.



ORDERING INFORMATION

Avantor ACE Method Development Kits are available in a range of dimensions and particle sizes (1.7 μm , 3 μm or 5 μm) to fit UHPLC and traditional HPLC systems.

Have questions? Reach out to our chromatography technical support team at chromsupport@avantorsciences.com or by calling **+44 (0) 118 930 3660**. Visit avantorsciences.com or ask your Avantor representative directly. We can help you with:

- Kit selection
- Product codes and dimensions
- Availability of guard columns and matching columns for scale-up

READY TO STREAMLINE YOUR NEXT METHOD?

Just take these steps:

- Begin with a structured screening strategy.
- Let selectivity lead the way.
- Use tools designed by chromatographers, for chromatographers.

Avantor ACE Method Development Kits deliver speed, reproducibility and confidence so you can focus on results, not revisions.



avantorsciences.com

Avantor® is a leading global provider of mission critical products and services to customers in the biopharma, healthcare, education & government, and advanced technologies & applied materials industries. We operate in more than 30 countries and deliver an extensive portfolio of products and services. We set science in motion to create a better world. Trademarks are owned by Avantor Inc. unless otherwise noted. © 2026 Avantor Inc.